



Sense

Using the latest in pattern recognition technology, Sense empowers users with intuitive control of their prostheses.

FlexCell

The first flexible battery for myoelectric prostheses. Place batteries virtually anywhere in the socket with ease to craft better-looking and long-lasting prostheses.

Thanks for buying one of our products.

Here are the other devices we offer.
There's more on the back, too!
Hope to hear from you again.

From your friends at



Check out our social media here!



(443) 451-7175


sales@i-biomed.com

www.i-biomed.com




IBT Research Kit

A platform to develop and implement your ideal control strategy for myoelectric prostheses.



SINCE 2011
OVER 3000 PATIENTS SERVED!



Element Electrodes

The most compact suction electrodes in the industry. Simplify the fitting process with wireless gain adjustment.





V2P

Easily adjust between light and heavy pinch forces with a built-in slider.



Equilux

The first body-powered terminal device that allows users to switch between voluntary opening and voluntary closing modes with the flip of a hinge.



Retro Classic

A no-nonsense hook with ample grip force for the most rugged jobs.



ITAL, HSC

A ready-to-use universal socket that's tough enough for everyday use.

Self-suspension is achieved using an innovative Humeral Suspension Cuff (HSC) that can be easily adjusted for comfort or to allow for limb volume changes. The HSC can be ordered separately to be used to suspend a custom fabricated prosthesis.

Partner Products

We believe that when patients have access to more products that share our vision of innovation and great engineering, everyone benefits. Thanks to your support, our distribution network allows companies who partner with us to reach more patients in need.

From our friends at



Z0000084rev3

Point Digits by

Made from 3D printed titanium, these mechanical fingers feature the most anatomically accurate center-of-rotation. An innovative ratcheting mechanism allows for truly hands-free operation while supporting up to 150 lbs of force. Point Designs offers the Point Digit, Point Digit mini, Point Partial, and Point Thumb.



Locking Wrists

Fully adjustable wrists that can lock in a vast number of positions.



Custom Attachments

Modified tools to fit the 1/2"-20" thread that is standard on all wrists.

